

Percussion Massage Belt



User's Manual

Thank you for purchasing the Percussion Massage Belt by Wagan Tech. With normal care and proper treatment it will provide years of reliable service. Please read all operating instructions carefully before use.

FEATURES

- 2 Rhythmic Massagers
- 10 Intensity levels
- 20 Different Combinations
- Versatile: Use for Neck/Shoulders, Back, Lumbar, Legs and Posterior

WARNING

To reduce the risk of burns, fire, electric shock or injury to user, please follow the precautions below

- Always unplug the unit when not in use.
- Do not stand or sit on the unit
- Only use the unit for the purposes described in this operating manual.
- Keep power cords away from heated surfaces. Do not carry the unit by the power cord. Do not allow pins or moisture to get inside the plug.
- Do not use the unit in the bathroom or other damp, humid places, for this may cause electric shock or cause the unit to malfunction. Do not spill liquid onto the control panel surface.
- Do not use the unit on top of heated appliances, including electric blankets.

- Never operate the unit with a damaged power cord or plug. When unplugging the unit, remember to turn off the unit first. Do not use the unit with a transformer as it may lead to malfunction or electric shock.
- To avoid electric shock or injury, do not attempt to open or disassemble any part of the unit.
- Do not attempt to repair the power cord if it is damaged. Please contact our company for assistance as needed.

SAFETY PRECAUTIONS TO AVOID RISK OF INJURY

- Do not use the unit against bare skin: while thin clothing may increase effectiveness, exposing the skin directly to the massage may irritate the skin.
- Person with any of the conditions listed below or who are undergoing medical treatment should consult with their physicians before using the unit:
 - Pregnant women, people who are ill or in poor physical shape.
 - People who are prohibited from receiving a massage by a physician due to thrombosis, aneurism, acute varicose veins or other circulatory disorders.
 - People with irregular curvatures of the spine.
- While using the massager, if you start feeling sick or if the massage seems painful, please stop use it immediately.
- Do not fall asleep while using the unit. Do not use the unit after consuming alcohol.
- Do not allow children to play with the unit.
- Do not unplug the unit when it is still massaging.
- Unplug the unit immediately if there is a power outage.

UNIT INSTRUCTIONS:

Place the Percussion Massage Belt on your shoulders, back or legs as shown in the Massage Areas section.


1. Connect the cable to a standard 110V AC power supply. After a beep, the appliance will go into standby mode and the POWER LED will blink red.
2. Press the **(M)** (mode) button and the unit will switch to automatic mode. The "AUTO" indicator will light green. The massager will automatically cycle through all 20 different massage rhythm combinations over a 15 minute program.
3. Press the **(M)** button again and the unit will switch to manual mode. The "MANUAL" indicator will light red. The massager will run a single massage combination continuously for 15 minutes. You can select a different massage combination by pressing the **(M)** button repeatedly until getting the desired pattern.
4. When the massager is operating: Press the **(+)** and **(-)** buttons to adjust intensity levels. There are 10 levels for adjustment.
5. To stop the unit, press and hold the **(M)** button for 3 seconds, and the unit will turn off.

Note: For both Automatic and Manual modes the massager will run for 15 minutes, after which the timer will automatically shut it off. Allow the unit to cool for at least 15 minutes before reusing.


Note: The Percussion Massage Belt has a safety feature causing the massage motors to shut off if overheated. If one or both massage motors does not activate, allow unit to cool significantly before attempting to use again.

MASSAGE AREAS


The unit massages the shoulder, back, lumbar, legs or posterior with rhythmic percussion. See illustrations below for usage examples.




Neck/Shoulders




Upper Back



Lumbar (Lower Back)



Legs



Posterior

CLEANING AND MAINTENANCE

- Always unplug the unit before cleaning it. Never touch the power plug with wet hands.
- Wipe the unit with a soft, dry cloth. Do not use cloths containing any kind of chemical.
- Never submerge any part of the unit in liquid.
- If the unit is particularly dirty:
 1. Dampen a cloth in water or a 3–5% mild detergent solution.
 2. Wipe the dirty areas with the wet cloth.
 3. Clean the cloth in water and wring it out, then wipe away any remaining detergent off the surface of the unit.
 4. Wipe the areas with a soft, dry cloth to dry.
 5. Wait for the unit to become fully dry before use.
- Do not press the unit against plastic for a long time, since it may lead to discoloration.
- Please do not use alcohol or any solvents to clean the unit.

SPECIFICATIONS

Rated Voltage	110V
Frequency	60Hz
Power	60W

WAGAN Corp. Limited Warranty

The WAGAN Corp. warranty is limited to products sold only in the United States.

All Wagan Tech products must be registered within (30) days of purchase to activate its warranty. To register your product, please visit <http://tinyurl.com/waganwarranty>. Be sure to keep the original receipt as it will be required when returning a product under the warranty.

Warranty Duration: This product is warranted to the original purchaser for a period of one (1) year from the original purchase date, to be free of defects in material and workmanship. WAGAN Corporation disclaims any liability for consequential damages. In no event will WAGAN Corporation be responsible for any amount of damages beyond the amount paid for the product at retail.

Warranty Performance: During the above one (1) year warranty period, a product with a defect will be replaced with a comparable model when the product is returned to WAGAN Corporation with an original store receipt. The replacement product will be in warranty for the balance of the original one (1) year warranty period.

To return a defective item, please contact WAGAN Corporation at (800) 231-5806 to obtain a Returned Merchandise Authorization number (RMA#), and return instructions. Each item returned will require a separate RMA#. After you have received the RMA# and the return instructions from WAGAN Corporation, please follow the instructions and send the item with PREPAID SHIPPING, along with all of the required documentation, a complete explanation of the problem, your name, address and daytime phone number. WAGAN Corporation will, at its option, replace or repair the defective part.

A Returned Merchandise Authorization number (RMA#) is REQUIRED when sending in any defective item. WAGAN Corporation is not responsible for any item(s) returned without an official Returned Merchandise Authorization number. The item(s) must be returned with prepaid shipping. WAGAN Corporation is not responsible for any shipping charges incurred in returning the item(s) back to the company for repair or replacement. This warranty is void if the product has been damaged by accident, in shipment, unreasonable use, misuse, neglect, improper service, commercial use, repairs by unauthorized personnel or other causes not arising out of defects in materials or workmanship. This warranty does not extend to any units which have been used in violation of written instructions furnished.

Warranty Disclaimers: This warranty is in lieu of all warranties expressed or implied and no representative or person is authorized to assume any other liability in connection with the sale of our products. There shall be no claims for defects or failure of performance or product failure under any theory of tort, contract or commercial law including, but not limited to negligence, gross negligence, strict liability, breach of warranty and breach of contract.

©2008 WAGAN Corp. REV2012



www.wagan.com
customerservice@wagan.com
USA & CAN Toll Free: 1.800.231.5806
Tel: +1.510.471.9221

31088 San Clemente Street
Hayward, CA 94544
U.S.A.



Intertek

4001459

©2011 Wagan Corporation. All Rights Reserved
Wagan Tech and wagan.com are trademarks used by Wagan Corporation
REV20110928